

Protocols for ET recombination (original protocol by F. Stewart at <http://www-db.embl-heidelberg.de/jss/servlet/de.embl.bk.wwwTools.GroupLeftEMBL/ExternalInfo/stewart/ETprotocols.html>).

(Modified by Alex Nechiporuk, March 2009)

Oligo design

* The 5' end (the homology arm) - choose 42 or more (I usually choose about 50) nts for the homology arms from the target DNA sequence simply according to where you want to insert the PCR fragment. I just order my primers de-salted, but if you are extra paranoid get them purified. Original protocol requires primer purification, I skip that step altogether and go to PCR.

33 ul dH₂O
5 ul 10 x PCR reaction buffer
5 ul 2.5 mM dNTP (10 x)
1.5 ul upper oligo
1.5 ul lower oligo
2 ul template
0.5 ul Taq polymerase (5 U/ul)

* Annealing temperature is important, usually 60°C - 62°C is optimal. Purify the PCR products by Qiagen column and elute with 2x 45 ul dH₂O.

* Add 10x Dpn 1 buffer and 2 ul Dpn 1 (NEB) and digest (to eliminate template DNA) for 1 hour.

* Extract with Phenol:CHCl₃ once and once with chloroform; add NaAc, pH=5.2 to a final concentration 0.3M and precipitate with ethanol; redissolve in 5 ul dH₂O (from 50 ul original PCR products, about 0.3 ug/ul).

Getting pBAD $\alpha\beta\gamma$ plasmid into your BAC clone

- Transform your BAC host with the ET expression plasmid (pBAD $\alpha\beta\gamma$) using electroporation procedure described below.
- Plate on 12.5 ug/ml chl and 100 ug/ml amp (if using pBAD $\alpha\beta\gamma$ that is amp resistant).
- Pick 3-6 single colonies and grow in 5 ml LB medium overnight.
- Extract using alkaline lyses and run on a gel against pBAD $\alpha\beta\gamma$ plasmid to check for plasmid integrity.

Making electrocompetent cells from your BAC culture (with or without pBAD $\alpha\beta\gamma$)

* Transfer 0.7ml into 70 ml of LB medium (without glucose!) and grow them at 37°C with shaking.

* Prepare 10% glycerol with dH₂O, and cool down on ice for at least 3 hours before using.

* When the cells reach OD₆₀₀ = 0.1-0.15, add 0.7ml 10% L-arabinose to induce ET protein expression (need to do this step only when BAC contains pBAD $\alpha\beta\gamma$).

* After a further 45-60 minutes, the cells should be at OD₆₀₀ of 0.3-0.4. Harvest.

* make sure the centrifuge and SA600 or SS34 (Sorvall) rotor is very cold by centrifuging for 10 min, 0°C at 4,000 rpm (I use regular 50 ml conical tubes with 600TC Sorvall rotor; Beckman has an analogous rotor for 50 ml conical tubes).

* Spin 35 mls cells for 10 min at 7,000 rpm at 0°C. Put the other 35 mls on ice.

* Pour away the supernatant, add the second 35 mls and respin.

* Pour away supernatant, put tube on ice, resuspend cells in 5 ml ice cold 10% glycerol with an ice cold 5ml pipette. Add a further 25 mls and centrifuge.

* Repeat the above step twice.

* Pour away supernatant and immediately dry the tube out with Kleenex tissue taking care not to touch the pellet.

* Resuspend the cells in the remaining liquid (you should have a little more than 100ul final resuspended volume).

* Transfer 50 ul of cells into each pre-cooled eppendorf tube and freeze in liquid N₂ or use immediately.

All the steps should be done as cold as possible, always on ice! The glass pipettes should be cooled once or twice by pipetting cold 10% glycerol up and down before pipetting cells (I found that glass pipettes work much better than plastic).

Transformation (same conditions for both pBAD $\alpha\beta\gamma$ and PCR product)

* (If cuvettes were reused) Wash cuvettes at least 10 times with dH₂O, precool them on ice for least 5 min.

* Thaw competent cells on ice and add 1 ul of PCR product (should be at least 0.1ug, 0.3 ug works better).

* Electroporate the cells at 2.3 kV (Bio-Rad Gene Pulser), 25 uF with Pulse controller set to 200 ohms)

* Add 1 ml of LB medium and transfer back into the eppendorf tube.

* Incubate at 37°C for 1 to 1.5 hours with shaking.

* Plate 100ul on chl/amp/kan (12.5/100/12.5 ug/ml) plates, spin the rest, resuspend in 100ul and plate that.

Genotyping

- Colonies may have an uneven appearance and that is usually a good indicator that reaction worked.
- Pick a single colony into 50 ul of LB containing chl/kan and use 2 ul of that in 20 ul PCR genotyping (no need to grow). I usually genotype between 24 and 48 colonies.
- Store the cultures at 4°C.
- Once a positive clone is identified (I've had a variable luck with >95% of all clones being positive or as little as 5%), pick into 100 ml of LB with chl/kan for midiprep. Without the amp bugs will rapidly lose pBAD $\alpha\beta\gamma$ plasmid.

BAC midiprep

- I use Clontech AX 100 columns, they provide a protocol for extracting low-copy plasmids. The typical yield is 25-50 ug from 100 ml of culture.
- At the end, resuspend the pellet in 50 ul of TE⁻⁴ (Tris HCl, pH=8.0 and **0.1** mM EDTA).

- Desalt the prep as follows: float Millipore 0.025 um membrane (cat#VSWP04700) shiny side up in 60 mm Petri dish filled about half-way with TE⁻⁴. Carefully pipet BAC prep on top of the membrane (will stay in a single drop), cover and let float for at least 2 hours. Collect with a pipet.
- Inject 20-120 ng/ul of supercoiled BAC diluted in water to check for transient expression.
- Rip the rewards...

References:

Zhang Y, Buchholz F, Muyrers JP, Stewart AF. A new logic for DNA engineering using recombination in Escherichia coli. Nat Genet. 1998 Oct;20(2):123-8.

Muyrers JP, Zhang Y, Testa G, Stewart AF. Rapid modification of bacterial artificial chromosomes by ET-recombination. Nucleic Acids Res. 1999 Mar 15;27(6):1555-7.